



SURGIPLUS

LOW FOAMING MULTI ENZYME CLEANER

Powerful pre soak surgical instrument and scope cleaner.

- Low-foaming solution designed for automated & Bath cleaning systems.
- Unparalleled combination of speed and efficacy
- Proprietary combination enzyme concentrate.
- Fully effective against proteins, carbohydrates, lipids and mucopolysaccharides.
- Unique stable formulation.
- Designed for use during ultrasonic or manual cleaning.
- Compatible with surgical, medical and dental instruments including flexible and rigid endoscopes.

Surgiplus is a unique combination of enzymes which effectively remove blood, tissue, mucous and other protein rich body fluids, from stainless steel, tungsten carbide, aluminium, brass, glass, plastic materials, and delicate micro surgical instruments, flexible and rigid scopes.

Surgiplus is use as per soak surgical instrument and scope cleaner.

Surgiplus immediately works to break down blood, fat and protein and saves reprocessing time by initiating the pre-cleaning process. For most application there is no need to add further cleaning ingredients for pre-soaking or manual cleaning.

Surgiplus save time, lowers inventory level, cuts the cost of consumables and improves turn around time.

Advantages :

- Provides quick and trough cleaning.
- Facilitates effective cleaning and rapid turn around.
- Immediately begins to break down blood, tissue, mucous and other protein rick body fluids.
- No other cleaning ingredients required.

Recommended Uses :

- Medical / surgical instruments.
- Respiratory therapy equipment.
- Evacuation system.
- Stainless steel, tungsten carbide, aluminium, brass, plastic materials.
- Delicate micro surgical instruments.
- Flexible and rigid scopes.

Direction For Use :

- Light biological contamination : 5 ml / ltr. of water-Soak time 2 to 5 min.
- Moderate biological contamination : 7.5 ml / ltr. of water-Soak time 8 to 10 min.
- Heavy biological contamination : 10 ml / ltr. of water-Soak time 10 min.
- Dired biological contamination : 20 ml / ltr. of water-Soak time 20 min.

Specifications :

Appearance	: Light Yellow Liquid.	Solubility	: Complete with water.
Rinsing property	: Excellent.	Foaming	: Low foam.
pH	: Neutral.	Biodegradable	: 100% Biodegradable.
Detergency	: Excellent.	Corrosion Factor	: Non-corrosive.

YOUR SURGICAL INSTRUMENT WILL LOOK NEWER AND PERFORM LONGER



Rudraksh Inc.

Mfgs. of Surgical Chemicals

13, Daji Colony, Nr. Vijaynagar Petrol Pump, VijayNagar, Naranpura, Ahmedabad - 380 013.
Mobile : +91 98985 57359 E-mail : rudrakshinc@yahoo.co.in